



### CONSTRUCTION OF STANDARD CURVES using GBlock

1. Create an account in [www.idtdna.com](http://www.idtdna.com)
2. Go to products and services
3. Select gBlocks® Gene Fragments. The price depends on the size of the DNA sequence
4. Select Order
5. Insert name of the sequence and the length of the synthetic DNA you want.

For HAdV

**INTEGRATED DNA TECHNOLOGIES**  
SPECIFICATION SHEET WWW.IDTDNA.COM

13-Nov-2013 Order No. **2274389**  
Ref. No. **66046805**

**Name - Standard HAdV** gBlocks® Gene Fragments 273 base pairs

5' - ATG ATG CCG CAA TGG TCT TAC ATG CAC ATC GCC GGG CAG GAC GCC TCG GAG TAT CTG AGC CCG GGC  
 ACA CAC ACA CAC CTG GTG CAA TTT GCC CGC GCC ACC GAT ACG TAC TTC AGC CTG GGG AAC AAG TTC AGA  
 AAT CCC GCT GCG ATT CGT GCC AGT CGA CCG CGA GGA CAC CGC TTA TTC TTA CAA AGT GCG CTT TAC GCT  
 GGC CGT GGG CGA CAA CCG GGT GTT GGA CAT GGC CAG CAC CTA CTT TGA CAT CCG CGG CGT GCT GGA TCG  
 - 3'

For MS2

**INTEGRATED DNA TECHNOLOGIES**  
SPECIFICATION SHEET WWW.IDTDNA.COM

01-Oct-2015 Order No. **2476419**  
Ref. No. **69797476**

**Name - MS2** gBlocks® Gene Fragments 346 base pairs

5' - CGT CGT AAG GTG CCT ACA AGC GAA GTG GGT CAT CGT GGG GTC GCC CGT ACG AGG AGA AAG CCG GTT  
 TCG GCT TCT CCC TCG ACG CAC GCT CCT GCT ACA GCC TCT TCC CTG TAA GCC AAA ACT TGA CTT ACA TCG  
 AAG TGC CGC AGA ACG TTG CGA ACC GGG CGT CGA CCG AAG TCC TGC AAA AGG TCA CCC AGG GTA ATT TTA  
 ACC TTG GTG TTG CTT TAG CAG AGG CCA GGT CGA CAG CCT CAC AAC TCG CGA CGC AAA CCA TTG CGC TCG  
 TGA AGG CGT ACA CTG CCG CTC GTC GCG GTA ATT GGC GCC AGG CGC TCC GCT ACC TTG CCC TAA ACG AAC  
 TGT C - 3'

6. When the stock arrive follow instruction for resuspension
7. Centrifuge the tube for 3-5 sec at a minimum of 3000 x g to ensure the material is in the bottom of the tube.
8. Add TE to reach a final concentration of 10 ng/ul. In the sheet attach with the product you can see the amount of DNA in ng.



**Properties**

Length:	346
Amount Delivered:	500ng
GC Content:	58.09%
Molecular Weight:	213696.6
fmoles/ng:	4.68
µg/OD <sub>260</sub> :	50

9. Vortex briefly.
10. Incubate at 50° C for 20 minutes.
11. Briefly vortex and centrifuge.
12. The tube resuspended have to be at -20°C
13. Quantify the gblock to know the exact amount of DNA per ul. We recommend Qubit R to do that. Normally less than 10 ng/ul.
14. Insert the data in the following excel sheet (double click on the table)

<b>Concentrations</b>					
qubit measure	7,29	ng/ul			
molecular weight	213696,6				
length in base pairs	273				
<b>Standard calculation</b>					
length	273	pb			
plasmid concentration	7,29	ng/µl			
Molecular weight	213696,60	u.m.a			
concentration genomic copies /gr	2,82E+18	cg/gr			
ng for 1E+11 genomic copies	35,49	ng			
volume that has 1E+11 genomic copies	4,87	µl			
volumen of TE buffer	995,13	ul			
stock concentration / ml	1,00E+11	1,00E+10	1,00E+09	1,00E+08	1,00E+07
genomic copies in 10 ul of reaction	1,00E+09	1,00E+08	1,00E+07	1,00E+06	1,00E+05
genomic copies in 5 ul of reaction	5,00E+08	5,00E+07	5,00E+06	5,00E+05	5,00E+04

15. The volume that has 1E+E11 genomic copies is diluted with the complement of TE to get 1ml



16. Serially dilute the DNA in order to obtain dilutions per 10 or 5  $\mu$ l. See the table in the excel sheet.
17. Test each dilution by triplicate in a qPCR reaction. Make sure signal (fluorescence, Ct values < 40) appear after  $5 \times 10^0$  copies/5  $\mu$ l for RNA or  $1 \times 10^0$  copies/10  $\mu$ l for DNA triplicates. If the signal appears at lower dilutions, the standard curve is not well constructed.

## Annex

### Sequences of the GBlock for the qPCR of Human adenoviruses and the bacteriophage MS2

>HAdV

```
ATGATGCCGCAATGGTCTTACATGCACATCGCCGGGCAGGACGCCTCGGAGTATCTGAGC
CCGGGCACACACACACACCTGGTGAATTTGCCGCGCCACCGATACGTACTTCAGCCTG
GGGAACAAGTTTCAGAAATCCCGCTGCGATTCGTGCCAGTCGACCGCGAGGACACCGCTTA
TTCTTACAAAGTGCGCTTTACGCTGGCCGTGGGCGACAACCGGGTGTGGACATGGCCAG
CACCTACTTTGACATCCGCGGCGTGCTGGATCG
```

>MS2

```
CGTCGTAAGGTGCCTACAAGCGAAGTGGGTCATCGTGGGGTCGCCCCGTACGAGGAGAAA
GCCGGTTTCGGCTTCTCCCTCGACGCACGCTCCTGCTACAGCCTCTTCCCTGTAAGCCAA
AACTTGACTTACATCGAAGTGCCGCAGAACGTTGCGAACCAGGGCGTCGACCGAAGTCCTG
CAAAAGGTCACCCAGGGTAATTTTAACTTGGTGTGGCTTTAGCAGAGGCCAGGTCGACAG
CCTCACAACCTCGCGACGCAAACCAATTGCGCTCGTGAAGGCGTACACTGCCGCTCGTCGCG
GTAATTGGCGCCAGGCGCTCCGCTACCTTGCCCTAAACGAACTGTC
```