

## Quantitative PCR SOPs for quantification of viruses in irrigation water

This SOP describes the performance of specific qPCR for 6 DNA viruses (human adenovirus (HAdV), JC polyomavirus (JCPyV), porcine adenovirus (PAdV), bovine polyomavirus (BPyV), ovine polyomavirus (OPyV) and avian parvovirus (Ch/TyPV)) and 1 RNA virus used as process control, MS2 phage.

All qPCRs have been tested by using the ABI PRISM®HID 7700 SDS from Applied Biosystems and the Mx3000P from Stratagene. However, several equipments can be used for running this SOP. The reactions take place in a 96-well optical reaction plate (Applied Biosystems. Cat. No. 4306737 for 20 units or 43426659 for 500 units) covered with ABI Prism™ optical adhesive covers (Applied Biosystems. Cat. No. 4311971 for 100 units or 4360954 for 25 units) and the optical adhesive cover starter kit (Cat. No 4313663). Optical caps and its corresponding installing tool can also be used (Cat. No 4323032 and 4330015).

Standards curves used in qPCR and qRT-PCR are generated by using serial dilutions of known amounts of a synthetic DNA construct containing the target region of the PCR. We strongly recommend storing the reference suspensions into DNA low binding tubes (Eppendorf Cat. No. 0030108-035 for 0.5 ml and Cat. No. 0030108-051 for 1.5 ml).

### Reagents:

**For qPCR:** TaqMan Environmental PCR Master Mix (Applied Biosystems. Part nº 4396838) supplied in a 2x concentration.

**For qRT-PCR:** RNA Ultrasense™ One-Step Quantitative RT-PCR System (Invitrogen Cat.No.: 11732-927).

### Mastermix preparation:

Prepare the mix in a clean separated area following the instructions provided:

qPCR	DNA	RNA
Environmental Mix 2x	12.5 µl	-----
Ultrasense Mix 5x	-----	5 µl
Primer Forward	1 µl	1 µl
Primer Reverse	1 µl	1 µl
Probe	0.5 µl	0.25 µl
ROX (dil 1/10)	-----	0.5 µl
Enzyme Mix	-----	1.25 µl
H <sub>2</sub> O	-----	11 µl
TOTAL VOLUME OF MIX	15 µl	20 µl
Sample	10 µl	5 µl

The inclusion of 0.10x more mastermix reaction is recommended to avoid inaccuracies during pipetting. Once the mix has been prepared aliquot 15 µl of mastermix for DNA viruses and 20 µl of mastermix for RNA viruses. An example of a qPCR template at the end of the SOP is provided.

### Primers and probes for qPCR and RT-qPCR

Virus	Ref	Name	Sequence (5'→3')	Working solutions*	Final PCR conc
HAdV	1, 2	AdF	CWTACATGCACATCKCSGG	22.5 µl	0.9 µM
		AdR	CRCGGGCRAAYTGACCAG	22.5 µl	0.9 µM
		AdP1	6-FAM-CCGGGCTCAGGTACTCCGAGGCGTCCT-BHQ1	11.25 µl	0.225 µM
JCPyV	3	JE3F	ATGTTTGCCAGTGATGATGAAAA	10 µl	0.4 µM
		JE3R	GGAAAGTCTTTAGGGTCTTCTACCTTT	10 µl	0.4 µM
		JE3P	6-FAM-AGGATCCCAACACTCTACCCACCTAAAAAGA-BHQ1	6 µl	0.12 µM
PAdV	4	Q-PAdV-F	AACGGCCGCTACTGCAAG	22.5 µl	0.9 µM
		Q-PAdV-R	AGCAGCAGGCTCTTGAGG	22.5 µl	0.9 µM
		Q-PAdV-P	6-FAM-CACATCCAGGTGCCGC-BHQ1	11.25 µl	0.225 µM
BPyV	5	QB-F1-1	CTAGATCCTACCCTCAAGGGAAT	10 µl	0.4 µM
		QB-R1-1	TTACTTGGATCTGGACACCAAC	10 µl	0.4 µM
		QB-P1-2	6-FAM-GACAAAGATGGTGTGTATCCTGTTGA-BHQ1	6 µl	0.12 µM
OPyV	6	Ov_qFI	TTCTCACAGCTGCAGACATTG	22.5 µl	0.9 µM
		Ov_qRI	TCCAATCTGGGCATAAGATTG	22.5 µl	0.9 µM
		Ov_qPI	6-FAM-ATGATTACCAAGCCAGACAGTGGG-BQH1	11.25 µl	0.225 µM
Ch/TyPV	7	Q-PaV-F	AGTCCACGAGATTGGCAACA	7.5 µl	0.3 µM
		Q-PaV-R	GCAGGTAAAGATTTTCACG	22.5 µl	0.9 µM
		Q-PaV-P	6-FAM-AATTATTCGAGATGGCGCCACG-BHQ1	12.5 µl	0.25 µM
MS2	8	pecson-2F	AAGGTGCCTACAAGCGAAGT	25 µl	1 µM
		pecson-2R	TTCGTTTAGGGCAAGGTAGC	25 µl	1 µM
		PecP-2	6-FAM-ATCGTGGGGTCGCCCCGTACG-BHQ1	25 µl	0.25 µM
NoVGI	9,10,11	QNIF4	CGCTGGATGCGNTTCCAT	12.5 µl	0.5 µM
		NV1LCR	CCTTAGACGCCATCATCTTAC	22.5 µl	0.9 µM
		TM9	6-FAM-TGGACAGGAGATCGC-MGB	25 µl	0.25 µM
NoVGII	12,13	QNIF2d	ATGTTTCAGRTGGATGAGRTTCTCWGA	12.5 µl	0.5 µM
		COG2R	TCGACGCCATCTTCATTCACA	22.5 µl	0.9 µM
		QNIFS	6-FAM-AGCACGTGGGAGGGCGATCG-TAMRA	25 µl	0.25 µM
HAV	14	HAV68	TCACCGCCGTTTGCCTAG	12.5 µl	0.5 µM
		HAV240	GGAGAGCCCTGGAAGAAAG	22.5 µl	0.9 µM
		HAV150	6-FAM-CCTGAACCTGCAGGAATTAA-MGB	25 µl	0.25 µM
HEV	15	JVHEVF	GGTGGTTTCTGGGGTGAC	25 µl	0.1 µM
		JVHEVR	AGGGGTTGGTTGGATGAA	25 µl	0.1 µM
		JVHEVP	6-FAM-TGATTCTCAGCCCTTCGC-BHQ1	5 µl	0.25 µM
RoV	16	NSP3-F	ACCATCTWCACRTRACCCTCTATGAG	10 µl	0.4 µM
		NSP3-R	GGTCACATAACGCCCTATAGC	10 µl	0.4 µM
		NSP3-P	6-VIC-AGTAAAAAGCTAACACTGTCAAA-MGB	10 µl	0.4 µM

\* Volume of starting stock solution of primer or probe to prepare 100µl of working solution (adjust to 100ul with DNA-RNA free water). A starting stock solution of 100 µM for both primers and probe is assumed

### **Analysis of target viruses:**

Add 10 µl of the nucleic acid extraction for DNA viruses or 5 µl of the nucleic acid extraction for RNA viruses per well in a separate area. Run direct and a ten-fold dilution in purified water of each extracted sample in duplicate.

### **Negative process control (NPC):**

Add 10 µl of the nucleic acid extraction for DNA viruses or 5 µl of the nucleic acid extraction for RNA viruses per well from the negative process control extraction. If positive results are obtained samples belonging to that concentration batch should be not be used for calculations.

### **NTC:**

Two wells for non-template controls (NTCs) should be included. NTC prove that mix does not produce fluorescence. Add 10 µl of molecular grade water for DNA viruses and 5µl for RNA viruses. Non-template controls shouldn't present signal (Cq). If an amplification signal (distinguishable from background fluorescence) is present, this would indicate cross-contamination or nonspecific amplification.

### **Negative extraction control:**

Two wells containing the negative extraction controls should be included for each batch of extraction performed.

Both NTCs and negative extraction controls should be negative. If positive results are obtained in these controls all positive samples should be retested.

To avoid any further contamination with standard suspension is strongly recommended to add samples and negative controls firstly, cover this part of the microplate, add standard positive control in a separate area and cover this part of the plate . One option is to cut the adhesive cover in two parts, one for covering the samples, the other for covering the wells containing the standard suspension.

### **Standard suspensions (dsDNA):**

In a separated area add  $10^1, 10^2, 10^3, 10^4$  and  $10^5$  GC/10 µl per reaction for DNA viruses or  $5 \times 10^1, 5 \times 10^2, 5 \times 10^3, 5 \times 10^4$  and  $5 \times 10^5$  GC/5 µl for RNA viruses of the dsDNA control materials previously diluted in a suitable buffer (e.g. TE buffer). Use a micropipette exclusively used for this purpose. Add the reference suspension in duplicate or triplicate.

### **Enzymatic inhibition control:** (suggested)

It is advisable to evaluate potential enzymatic inhibition due to inhibitors present in each of the different environmental matrices evaluated. To do so, in representative samples of each matrix analyzed, an extra inhibition control could be used.

For DNA viruses:

1 µl of reference suspension containing 10<sup>5</sup>GC/10 µl in 10 µl of water

1 µl of reference suspension containing 10<sup>5</sup>GC/10 µl in 10 µl of sample

1 µl of reference suspension containing 10<sup>5</sup>GC/10 µl in 10 µl of sample diluted 1:10

For RNA viruses:

1 µl of reference suspension containing 10<sup>5</sup>GC/10 µl in 5 µl of water

1 µl of reference suspension containing 10<sup>5</sup>GC/10 µl in 5 µl of sample

1 µl of reference suspension containing 10<sup>5</sup>GC/10 µl in 5 µl of sample diluted 1:10

### **Process control:**

Extract 140 µl of the stock suspension of MS2 phage used to spike all processed samples. A quantification of MS2 per sample will be performed to estimate virus recovery/extraction efficiency. If negative results are obtained or concentration efficiency is lower than <0.1% results obtained for this sample should be discarded.

### **TaqMan Assay parameters:**

Perform the QPCR into an adequate system selecting the appropriate parameters:

Virus	Genome	qPCR cycles
HAdV	DNA	10min at 95°C, 40 cycles (15s at 95°C and 1min at 60°C)
JCPyV	DNA	10min at 95°C, 40 cycles (15s at 95°C and 1min at 60°C)
PAdV	DNA	10min at 95°C, 45 cycles (15s at 95°C, 20s at 55°C and 20s at 60°C)
BPyV	DNA	10min at 95°C, 45 cycles (15s at 95°C, 30s at 60°C)
OPyV	DNA	2min at 50°C, 10min at 95°C, 40 cycles (15s at 95°C, 1min at 60°C)

<b>Ch/TyPV</b>	DNA	2min at 50°C, 10min at 95°C, 40 cycles (15s at 95°C, 1min at 60°C)
<b>MS2</b>	RNA	Following the RT-PCR (1h at 55°C), 5min at 95°C, 40 cycles (15s at 95°C, 1min at 60°C and 1min at 65°C)
<b>NoVGI</b>	RNA	Following the RT-PCR (1h at 55°C), 5 min at 95°C, 45 cycles (15s at 95°C, 1 min at 60°C, 1min at 65°C)
<b>NOVGII</b>	RNA	Following the RT-PCR (1h at 55°C), 5 min at 95°C, 45 cycles (15s at 95°C, 1 min at 60°C, 1min at 65°C)
<b>HAV</b>	RNA	Following the RT-PCR (1h at 55°C), 5 min at 95°C, 45 cycles (15s at 95°C, 1 min at 60°C, 1min at 65°C)
<b>HEV</b>	RNA	Following the RT-PCR (30min at 50°C), 15 min at 95°C, 45 cycles (10s at 95°C, 20s at 55°C, 15s at 72°C)
<b>RoV</b>	RNA	Following the RT-PCR (30 min at 48°C), 10 min at 95°C, 40 cycles (15s at 95°C, 1 min at 60°C)
<b>MNV</b>	RNA	

Once the run is completed store data and obtain the results as described in the user's manual of the equipment used.

### **Quality control**

Standard precautions should be applied in all the manipulations in order to reduce the probability of sample contamination by amplified DNA molecules. Decontaminate the pipettes after each use. Use UV or cleaning products that can be obtained from your local pipette suppliers.

### **Analysis of results:**

The amount of DNA will be defined as the median of the data obtained after correcting the dilution factor when needed.

Take into account that DNA suspensions containing less than 1 GC/reaction shouldn't present C<sub>q</sub> values. Decimal dilutions of standard DNA or RNA suspensions should present C<sub>q</sub> values within approximately 3 C<sub>q</sub>s of difference between one dilution and the next ten-fold dilution. However, this may not be the case in environmental samples where direct dilutions sometimes show enzymatic inhibition (higher C<sub>q</sub> values than expected considering the 1/10 dilution results); these results presenting inhibition should not be considered for quantification, considering only the assays with 1/10 dilutions.

All amplification plots shall be checked to identify false-positive results (reactions with C<sub>q</sub> values not associated with exponential amplification) caused by high or uneven background signal. This shall be noted and results for any reactions affected in this way shall be regarded as negative.

The slope of the regression curve should be between 3.10 – 3.60 (corresponding to amplification efficiencies of ~90-110%). Points of a minimum of 3 dilutions of standard must be included in the regression curve. Curves with  $r^2$  values of  $<0.98$  or with slopes out of the range 3.10-3.60 should not be used for calculations.

### Interpretation of results

If higher inhibition is observed/ expected the analysis of a -2 dilution may be considered. In samples showing  $C_q > 38$ , if two replicates are positive, we may consider the data as a robust and confirmed result. However, if only one positive signal with  $C_q > 38$  is detected in one sample, the result will be considered as negative if no other assay to confirm the positive result is performed. It is advisable to repeat the assay or perform a nested-PCR assay on these samples to rule out non-specific amplification. If nested-PCR shows negative results the sample should be considered negative since no confirmation has been possible.

Examples:

D <sub>1</sub>	D <sub>2</sub>	Dil <sub>1</sub>	Dil <sub>2</sub>	Conclusion
35	35.4	38.2	38.7	ok, calculate mean value
35	<del>38</del>	38.2	38.7	probable mistake, calculate mean of the 3 valid data
<del>35</del>	<del>35.4</del>	35.2	36	potential inhibition, consider for quantification only 1/10 dilutions
39	38.6	no Ct	no Ct	ok, calculate mean value
38.2	no Ct	no Ct	no Ct	repeat, or confirm by nested-PCR
No Ct	no Ct	no Ct	39.3	repeat, or confirm by nested-PCR

D: direct, Dil: 1-fold dilution

Expression of results: The results will be expressed in the tables as the row Ct values obtained in each assay, with a blank in the wells not used for the calculations of the mean GC values. An additional column will show the mean value as GC/100 ml.

### References:

1. Hernroth BE, Conden-Hansson AC, Rehnstam-Holm AS, Girones R, Allard AK. Environmental factors influencing human viral pathogens and their potential indicator organisms in the blue mussel, *Mytilus edulis*: the first Scandinavian report. Appl. Environ. Microbiol. 2002. 68: 4523-33.
2. Bofill-Mas S, Albinana-Gimenez N, Clemente-Casares P, Hundesa A, Rodriguez-Manzano, J, Allard, A, Calvo, M, Girones R. Quantification and stability of human

- adenoviruses and polyomavirus JCPyV in wastewater matrices. *Appl. Environ. Microbiol.* 2006. 72: 7894-7896.
3. Pal A, Sirota L, Maudru T, Peden K, Lewis AM. Real-time PCR assays for the detection of virus-specific DNA in samples with mixed populations of polyomaviruses. *J. Virol. Methods.* 2006 Jul 135(1):32-42.
  4. Hundesa et al. Development of a qPCR assay for the quantification of porcine adenoviruses as an MST tool for swine fecal contamination in the environment. *J Virol Methods.* 2009. 158(1-2):130-135.
  5. Hundesa A, Bofill-Mas S, Maluquer de Motes C, Rodriguez-Manzano J, Bach A, Casas M, Girones R. Development of a quantitative PCR assay for the quantitation of bovine polyomavirus as a microbial source-tracking tool. *J Virol. Methods.* 2010. Feb 163(2):385-389.
  6. Rusiñol M., Carratalà A., Hundesa A., Bach A., Kern A., Vantarakis A., Girones R., Bofill-Mas S.. Description of a novel viral tool to identify and quantify ovine faecal pollution in the environment. *Sc Total Env* 2013. 458-460C:355-360.
  7. Carratalà A., Rusiñol M., Hundesa A., Biarnés M., Rodriguez-Manzano J., Vantarakis A., Kern A., Suñen E., Girones R., Bofill-Mas S. A novel tool to trace poultry fecal contamination in the environment: specific detection and quantification of chicken/turkey parvoviruses. *Appl. Environ. Microbiol.* 2012. 78(20):7496-9
  8. Brian M. Pecson, Luisa Valério Martin and Tamar Kohn. Quantitative PCR for determining the infectivity of bacteriophage MS2 upon inactivation by heat, UV-B Radiation, and singlet oxygen: advantages and limitations of an enzymatic treatment to reduce false-positive results.. *Appl. Environ. Microbiol.* Sept. 2009, p. 5544-5554-33.
  9. Ref de Da Silva A.K., Le Saux J.C., Parnaudeau S., Pommepuy M., Elimelech M., Le Guyader F.S. Evaluation of removal of noroviruses during wastewater treatment, using real-time reverse transcription-PCR: Different behaviors of genogroups I and II. *Appl. Environ. Microbiol.* 2007, 73 pp. 7891–7897
  10. Svraka S., Duizer E., Vennema H., De Bruin E., Van Der Veer B., Dorresteyn B. et al. Etiological role of viruses in outbreaks of acute gastroenteritis in The Netherlands from 1994 through 2005. *J. Clin. Microbiol.* 2007, 45 pp. 1389–1394
  11. Hoehne M., Schreier E. Detection of Norovirus genogroup I and II by multiplex real-time RT-PCR using a 3'-minor groove binder-DNA probe. *BMC. Infect. Dis.* 2006, 6:69
  12. Loisy F., Atmar R.L., Guillon P., Le Cann P., Pommepuy M., Le Guyader F.S. Real-time RT-PCR for norovirus screening in shellfish. *J. Virol. Methods.* 2005, 123 pp. 1–7

13. Kageyama T., Kojima S., Shinohara M., Uchida K., Fukushi S., Hoshino F.B. et al. Broadly reactive and highly sensitive assay for Norwalk-like viruses based on real-time quantitative reverse transcription-PCR. *J. Clin. Microbiol.* 2003, 41 pp. 1548–1557
14. Costafreda M.I., Bosch A., Pintó R.M. Development, evaluation, and standardization of a real-time TaqMan reverse transcription-PCR assay for quantification of hepatitis A virus in clinical and shellfish samples. *Appl. Environ. Microbiol.* 2006, 72 pp. 3846–3855.
15. Jothikumar N., Cromeans T.L, Robertson B. H, Meng X.J., Hill V.R. A broadly reactive one-step real-time RT-PCR assay for rapid and sensitive detection of hepatitis E virus. *Journal of Virological Methods.* 2006, 131 pp. 65-71.
16. Zeng SQ, Halkosalo A, Salminen M, Szakal ED, Puustinen L, Vesikari T. One-step quantitative RT-PCR for the detection of rotavirus in acute gastroenteritis. *J Virol Methods.* 2008 153, 238-40.



### Example plate layout (HAdV)

HAdV dsDNA 10 <sup>1</sup>	HAdV dsDNA 10 <sup>1</sup>	HAdV dsDNA 10 <sup>2</sup>	HAdV dsDNA 10 <sup>2</sup>	HAdV dsDNA 10 <sup>3</sup>	HAdV dsDNA 10 <sup>3</sup>	HAdV dsDNA 10 <sup>4</sup>	HAdV dsDNA 10 <sup>4</sup>	HAdV dsDNA 10 <sup>5</sup>	HAdV dsDNA 10 <sup>5</sup>		
Cut with the scissors the optical adhesive and cover the sample wells to avoid any further contamination											
										CNP	CNP
Test sample 1 (neat)	Test sample 1 (neat)	Test sample 1 (-1)	Test sample 1 (-1)	Test sample 2 (neat)	Test sample 2 (neat)	Test sample 2 (-1)	Test sample 2 (-1)	Test sample 3 (neat)	Test sample 3 (neat)	Test sample 3 (-1)	Test sample 3 (-1)
Test sample 4 (neat)	Test sample 4 (neat)	Test sample 4 (-1)	Test sample 4 (-1)	Test sample 5 (neat)	Test sample 5 (neat)	Test sample 5 (-1)	Test sample 5 (-1)	Test sample 6 (neat)	Test sample 6 (neat)	Test sample 6 (-1)	Test sample 6 (-1)
Test sample 7 (neat)	Test sample 7 (neat)	Test sample 7 (-1)	Test sample 7 (-1)	Test sample 8 (neat)	Test sample 8 (neat)	Test sample 8 (-1)	Test sample 8 (-1)	NTC	NTC	CNE	CNE

Standard suspension: 10µl for DNA viruses, 5µl for RNA viruses, 15 µl mastermix for DNA viruses and 20 µl for RNA viruses
Nucleic acid extraction: 10µl for DNA viruses, 5µl for RNA viruses, 15 µl mastermix for DNA viruses and 20 µl for RNA viruses
CNP: 10µl for DNA viruses, 5µl for RNA viruses, 15 µl mastermix for DNA viruses and 20 µl for RNA viruses using the control process NA extraction
NTC: 10µl of molecular grade water for DNA viruses, 5µl of molecular grade water for RNA viruses, 15 µl mastermix for DNA viruses and 20 µl for RNA viruses
CNE: 10µl of negative extraction control for DNA viruses, 5µl of negative extraction control for RNA viruses, 15 µl mastermix for DNA viruses and 20 µl for RNA viruses

Example plate layout for process control (suggested process control MS2 phage)

MS2 dsDNA 10 <sup>1</sup>	MS2 dsDNA 10 <sup>1</sup>	MS2 dsDNA 10 <sup>2</sup>	MS2 dsDNA 10 <sup>2</sup>	MS2 dsDNA 10 <sup>3</sup>	MS2 dsDNA 10 <sup>3</sup>	MS2 dsDNA 10 <sup>4</sup>	MS2 dsDNA 10 <sup>4</sup>	MS2 dsDNA 10 <sup>5</sup>	MS2 dsDNA 10 <sup>5</sup>		
Test sample 1 (neat)	Test sample 1 (-1)	Test sample 2 (neat)	Test sample 2 (-1)	Test sample 3 (neat)	Test sample 3 (-1)	Test sample 4 (neat)	Test sample 4 (-1)	Test sample 5 (neat)	Test sample 5 (-1)	Test sample 6 (neat)	Test sample 6 (-1)
Test sample 7 (neat)	Test sample 7 (-1)	Test sample 8 (neat)	Test sample 8 (-1)	NTC	NTC	CNE	CNE				

Standard suspension: 5 µl for dsDNA MS2 phage standard suspension and 20 µl of mastermix
Nucleic acid extraction: 5 µl for extracted RNA viruses and 20 µl of mastermix
NTC: 5 µl of molecular grade water for RNA viruses and 20 µl for RNA viruses
CNE: 5 µl of negative extraction control for RNA viruses and 20 µl for RNA viruses